


Specification

TE buffer (1X) pH 8.0 for molecular biology

A0386

Physical Description:	Liquid
Product Code:	A0386
Product Name:	TE buffer (1X) pH 8.0 for molecular biology
Specifications:	<p>DNases/RNases/Proteases: not detectable pH (20°C; adjusted with HCl): 8.0 ± 0.1</p> <p>Composition:</p> <p>EDTA · Na₂: 1 mM Tris: 10 mM</p>
WGK:	nwg
Storage:	RT
CS:	38221900
	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340