

Specification

Tween® 20 for molecular biology

A4974

Solubility:	100 g/L (H ₂ O)
Physical Description:	Liquid
Product Code:	A4974
Product Name:	Tween® 20 for molecular biology
Headline Comment:	®trademark of ICI America Inc.
Specifications:	DNases/RNases/Proteases: not detectable Density (d 20°C/4°C): 1.095 - 1.105 Hydroxyl-No.: approx. 100 Saponification No.: 40 - 50
WGK:	1
Storage:	RT
Molecular Formula:	C ₅₈ H ₁₁₄ O ₂₆
M:	1227.72 g/mol
CAS:	9005-64-5
EINECS:	500-018-3
CS:	34024200

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 07 228 16476 • Register court Darmstadt HRB Nr. 7340