

Specification

Urea for molecular biology

A1049

Solubility:	1000 g/L (H ₂ O)
Physical Description:	Solid
Product Code:	A1049
Product Name:	Urea for molecular biology
Specifications:	DNases/RNases/Proteases: not detectable Assay (from N): min. 99.0 % pH (8 M; H ₂ O; 20°C): 7.5 - 10.5 Heavy metals (as Pb): max. 0.001 % Chloride: max. 0.0005 % Sulfate: max. 0.001 % Fe: max. 0.0005 %
WGK:	1
Storage:	RT
Molecular Formula:	NH ₂ CONH ₂
M:	60.06 g/mol
CAS:	57-13-6
EINECS:	200-315-5
CS:	31021010

AppliChem GmbH