



PharmaTainer – Product Specification Sheet

10L PET Pharmatainer (UniPack)	Catalogue number 110000-01
---------------------------------------	-----------------------------------

Bottle part number	110000-single	USP Class VI / EP	ADCF	DMF No.
Bottle material	Cellon resin code 001	Yes	Yes	15205
Nominal volume	10,000ml	Tolerance +/-2.1mm	Accuracy +/- 2%	Shape (footprint) Square
Brim capacity	11,980ml			
Bottle dimensions w/o cap	H395 x W212 x D212mm			
Bottle weight w/o cap	556g			
Bottle thread O.D.	63.77mm	+/-0.3mm		
Bottle neck I.D.	54.5mm	+/-0.30mm		
Graduations	1000ml – 10,000ml	500ml intervals		Endotoxin
Sterility	C60 γ irradiated	SAL 10-6		No detectable Endotoxin Certified by Lot
Bottle decoration	Graduation scale	Lot number	Expiration date	Data matrix
Bottle colour	Clear			
Traceability	Lot number	Data matrix	Serial number	

Closure part number	300070C1-01	USP Class VI / EP	ADCF	DMF No.
Closure material	Cellon resin code 003	Yes	Yes	1646
		Tolerance +/-3%		
Closure weight	23.4g			
Closure dimensions	H 28.0mm O.D. 70.0mm	+/- 0.3mm		
Closure colour	Natural			

Packaging materials	Carton	Level 3 Carton liner	Level 2	Level 1 Sterile barrier
	4mm single corrugate, top and base protection sheet	50 μ LDPE monolayer film tape sealed bag	100 μ LDPE monolayer film Heat sealed bag	120 μ PE/PA/PE multilayer film vacuum bag heat sealed

Carton dimensions	Carton weight	Units per carton	Units per sterile barrier
L 458, W 248, H 458mm	2.13kg	2 units	1 (UniPack 1)

	Dimensions L x W x H	Weight (inc.pallet)	Cartons per pallet	Units per pallet
Euro pallet	1200x 800x1982mm	78.2kg	28	56
Industry pallet	1200x1000x1782mm	85kg	34	68

VERSION NO	003	28/02/2012
-------------------	------------	-------------------