








# CELLON

Revolutions in Cell & Tissue Culture

Bouchon Multiports Autoclavables pour Flacons PharmaTainer

Cap	Catalog Number	Description	External			Internal			PharmaTainer
			Nb of ports	Out. Diameter	In. Diameter	Nb of ports	Out. Diameter	In. Diameter	
	4048CP2-50	Cap 48 mm 2 ports	2	8,3 mm	6,4 mm	2	8,3 mm	6,1 mm	For PharmaTainer PET or PC from 125 ml to 2 liters.
	4048CP3-50	Cap 48 mm 3 ports	2	8,3 mm	6,4 mm	2	8,3 mm	6,1 mm	
			1	5,8 mm	4,3 mm	na	na	na	
	4148CP3-50	Cap 48 mm 3 ports	2	8,3 mm	6,4 mm	2	8,3 mm	6,1 mm	
			1	5,8 mm	4,3 mm	1	5,7 mm	3,6 mm	
	4070CP2-50	Cap 70 mm 2 ports	2	8,3 mm	6,4 mm	2	8,3 mm	6,1 mm	
	4070CP3-50	Cap 70 mm 3 ports	2	8,3 mm	6,4 mm	2	8,3 mm	6,1 mm	
			1	5,8 mm	4,3 mm	na	na	na	
	4170CP3-50	Cap 70 mm 3 ports	2	8,3 mm	6,4 mm	2	8,3 mm	6,1 mm	
			1	5,8 mm	4,3 mm	1	5,7 mm	3,6 mm	
	4070CP4-50	Cap 70 mm 4 ports	4	10,6 mm	8,3 mm	4	10,6 mm	8,0 mm	

