

## FICHE TECHNIQUE

---

<b>Référence</b>	351152
<b>Désignation</b>	Corning® FluoroBlok™ 24-Multiwell Insert Systems, PET Membrane



### Description :

- ✓ Plate Format : 24-well
- ✓ Cell Growth Area : 0,3 cm<sup>2</sup>
- ✓ Surface Treatment : TC-Treated
- ✓ Sterile : Yes
- ✓ Membrane Material : Polyethylene terephthalate (PET)
- ✓ Membrane Pore Size : 8 μm
- ✓ Nominal Pore Density : 6+/-2 x 10<sup>4</sup> pores per cm<sup>2</sup>
- ✓ Membrane Diameter : 6,4 mm
- ✓ Storage : Room temperature

### Unité de vente :

48 inserts

---