



Brand: Honeywell Fluka
Product: 35256

Sodium hydroxide solution 1 mol/l, volumetric solution, Reag. Ph. Eur.

Parameter	Specification
factor	0.999-1.001
BAM/EMPA ref.-material	effektiv
BAM/EMPA Certificate no.	effektiv
Traceable to NIST SRM	effektiv
Produced using quality grade "puriss. p.a." tested raw materials.	

Test method: Acidimetric titration related to hydrochloric acid measured against EMPA/BAM reference material (given above)
Product-No 71363: Sodium carbonate or
Product-No 93440: TRIS (Tris(hydroxymethyl)aminoethane).

The factor (titer) is measured at 20°C.

BAM: Federal Institute for Materials Research and Testing
Germany
EMPA: Swiss Federal Laboratories for Materials Testing and Research
NIST: National Institute of Standards and Technology
USA