

CORNING	Life Sciences	Document Number: LSR00137	Rev.: 5
Description (Class): Reference, Customer Technical Data Sheet 35x10 mm Easy Grip Dish			Page: 1 of 2

APPLICATIONS:

- Growth of cells.
- Unimpaired observation of dish contents.
- Ideal for microphotography.

FEATURES:

- Flat, distortion free optics.
- Lids designed for optimal gas exchange.
- Stacking rings allow for easier stacking and handling.
- Frosted rim improves the handling of the dish. This feature facilitates the ability to pick up the dish without accidentally removing the lid to guarantee aseptic manipulation.
- Reliable vacuum-gas plasma treatment cleans and permanently and consistently modifies the cell-growth surface.
- Non-pyrogenic
- Packaged in peel-open medical-style bags.
- Durable construction for stable dish manipulation

REGULATORY COMPLIANCE:

Falcon® products are manufactured under the current ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

STERILITY:

Product is gamma irradiated

RESTRICTIONS:

- Refer to "Thermoplastics Properties Chart," available at:
<https://www.corning.com/catalog/cls/documents/compatibility-guides/CLS-DL-AN-226.pdf>


Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.

Description (Class):
Reference, Customer Technical Data Sheet
35x10 mm Easy Grip Dish



Page:
2 of 2

All dimensions in millimeters (mm)


- Assembled Dish and Lid height: 12.8 mm
- Effective growth area = 9.6 cm²
- Working Volume = 2.5 - 3.0 mL
- Actual dimensions: 40.28 mm O.D. x 6.17 mm
- Actual growth area: 11.78 cm²



LID



DISH



MATERIAL: CRYSTAL-GRADE POLYSTYRENE

ITEM	CAT. NO.	DESCRIPTION	DISH AND LID KEY DIMENSIONS											
			A	B	C	D	E	F	G	H	J	K	L	M
1	353001	35x 10 mm Easy Grip Dish. TC-treated. 20/sleeve, 500/case	37.13	35.61	35.00	37.90	39.14	11.05	9.56	1.01	40.28	38.78	6.17	.889
2	351008	35x 10 mm Easy Grip Dish Non-TC-treated. 20/sleeve, 200/case												
2	353801	35x 10 mm Easy Grip Dish. Corning® Primaria™ -reated. 20/sleeve, 500/case												

Warranty/Disclaimer: Unless otherwise specified, all products are For Research Use Only. Not intended for use in diagnostic or therapeutic procedures.