


FICHE TECHNIQUE

Référence	353227
Désignation	Falcon® 96-well Polystyrene Microplates



Description :

- ✓ Plate Format : 96-well
- ✓ Plate Color : Clear
- ✓ Well Bottom : Round
- ✓ Well Volume : 0,32 mL
- ✓ Cell Growth Area : 0,36 cm²
- ✓ Recommended Working Volume : 50-250 µL
- ✓ Surface Treatment : TC-Treated
- ✓ Sterile : Yes
- ✓ Lids Included : Oui
- ✓ Storage : Ambient storage conditions

Unité de vente :

50 Plates
