

Brand: Honeywell Fluka
Product: 35354

Sulfuric acid 0,5 mol/l, volumetric solution Reag. Ph. Eur.

Parameter	Specification
factor	0.999-1.001
BAM/EMPA ref.-material	effektiv
BAM/EMPA Certificate no.	effektiv
Traceable to NIST SRM	effektiv
Produced using quality grade "puriss. p.a." tested raw materials.	

Test method: Acidimetric titration related to
EMPA/BAM referencematerial (given above)
Product-No 71363: Sodium carbonate or
Product-No 93440: TRIS (Tris(hydroxymethyl)aminoethane).

The factor (titer) was measured at 20°C.

BAM: Federal Institute for Materials Research and Testing
Germany

EMPA: Swiss Federal Laboratories for Materials Testing and Research

NIST: National Institute of Standards and Technology
USA