

FICHE TECHNIQUE

Référence	353938
Désignation	Plaque insert

Description :

- ✓ Storage: Store at room temperature, away from direct sunlight.
- ✓ Automation compatible design
- ✓ Format compatible with most robots and fluid handling instruments
- ✓ Complete sample recovery
- ✓ The Falcon 96 Square Well, Angled-Bottom Plate features an angled bottom for more complete sample utilization
- ✓ Excellent reproducibility
- ✓ One-piece feeder tray enhances consistency in well-to-well monolayer growth
- ✓ Total assay flexibility—ideal for transport studies
- ✓ System can be used with many cell lines including Caco-2, MDCK, and LLC-PK1, for basal to apical or apical to basal measurements of drug transport
- ✓ Falcon® 96 Square Well Angled Bottom Insert Plate, with Lid, Sterile

Unité de vente :

5 unités
