


FICHE TECHNIQUE

Référence	353988
Désignation	Falcon® 384-well White Flat Bottom



Description :

- ✓ Falcon® 384-well White Flat Bottom TC-treated Microtest Microplate, with Lid, Sterile
- ✓ Brand Falcon®
- ✓ Plate Format 384-well
- ✓ Plate Color Clear
- ✓ Well Bottom Flat
- ✓ Well Volume 130 µL
- ✓ Cell Growth Area 10 mm²
- ✓ Recommended Working Volume 15-100 µL
- ✓ Surface Treatment TC-Treated
- ✓ Sterile Yes
- ✓ Lids Included Oui
- ✓ Storage Ambient storage conditions

Unité de vente :

50 unités
