

## FICHE TECHNIQUE

---

<b>Référence</b>	354259
<b>Désignation</b>	Corning® Laminin/Entactin High Concentration, 10.5 mg



### Description :

- ✓ Purity > 95% by SDS-PAGE
- ✓ Tested for ability to initiate neurite outgrowth of NG-108 Mouse neuroblastoma cells
- ✓ Found negative for bacteria, fungi, and mycoplasma
- ✓ Stable for at least three months at -70 °C. Avoid multiple freeze-thaws
- ✓ Brand Corning®
- ✓ Size 10.5 mg
- ✓ Source Engelbreth-Holm-Swarm lathrytic mouse tumor
- ✓ Formulation High concentration laminin/entactin complex frozen in Dulbecco's PBS

### Unité de vente :

1 unité

---