

FICHE TECHNIQUE

Référence	354265
Désignation	Corning® Collagen I



Description :

- ✓ Electrophoretically homogeneous (SDS-PAGE)
- ✓ Shown to promote attachment and spreading of HT-1080 human fibrosarcoma cells
- ✓ Tested and found negative for bacteria, fungi, and mycoplasma
- ✓ Brand : Corning®
- ✓ Size : 10 mg
- ✓ Source : Human placenta
- ✓ Formulation : Frozen in 2mM HCl

Unité de vente :

1 unité
