


## FICHE TECHNIQUE

---

<b>Référence</b>	354410
<b>Désignation</b>	Plaque 96 puits



---

### Description :

- ✓ Corning® BioCoat® Laminin 96-well Clear Flat Bottom TC-treated Microplate with Lid
- ✓ Brand Corning®
- ✓ Plate Format 96-well
- ✓ Plate Color Clear
- ✓ Well Bottom Flat
- ✓ Well Volume 370  $\mu$ L
- ✓ Cell Growth Area 0,32 cm<sup>2</sup>
- ✓ Sterile No
- ✓ Lids Included Oui
- ✓ Storage 2°C- 8°C
- ✓ Surface Coating Laminin

### Unité de vente :

5 unités

---