


## FICHE TECHNIQUE

---

<b>Référence</b>	354412
<b>Désignation</b>	Plaque culture 24 puits



---

### Description :

- ✓ Corning® BioCoat™ Laminin 24-well Clear Multiwell Plate is comprised of clear TC-treated polystyrene with a uniform application of mouse laminin
- ✓ Brand : Corning®
- ✓ Cell Growth Area : 2 cm<sup>2</sup>
- ✓ Storage : 2°- 8°C
- ✓ Surface Coating : Laminin

### Unité de vente :

5 unités

---