

## FICHE TECHNIQUE

---

<b>Référence</b>	354428
<b>Désignation</b>	Plaque culture 6 puits



---

### Description :

- ✓ Corning® BioCoat™ Collagen IV 6-well Clear Multiwell Plate is comprised of clear TC-treated polystyrene with a uniform application of mouse collagen type IV. Flat bottom
- ✓ Brand Corning®
- ✓ Storage 2° - 8°C
- ✓ Surface Coating Collagen IV

### Unité de vente :

5 unités

---