

## FICHE TECHNIQUE

---

<b>Référence</b>	354429
<b>Désignation</b>	Plaque 96 puits



---

### Description :

- ✓ Corning® BioCoat® Collagen IV 96-well Clear Flat Bottom TC-treated Microplate, with Lid
- ✓ Brand Corning®
- ✓ Plate Format 96-well
- ✓ Plate Color Clear
- ✓ Well Bottom Flat
- ✓ Well Volume 370  $\mu$ L
- ✓ Cell Growth Area 0,32 cm<sup>2</sup>
- ✓ Sterile No
- ✓ Lids Included Oui
- ✓ Storage 2°C- 8°C
- ✓ Surface Coating Collagen IV

### Unité de vente :

5 unités

---