

FICHE TECHNIQUE

Référence	354444
Désignation	Insert



Description :

- ✓ Brand Corning®
- ✓ Cell Growth Area 0.3 cm²
- ✓ Membrane Material Polyethylene terephthalate (PET)
- ✓ Membrane Pore Size 0.4 µm
- ✓ Membrane Diameter 6.4 mm
- ✓ Packaging Format 12 Inserts/Plate/Pk, 2 Plates/Cs
- ✓ Application Air-Liquid Interface (ALI) Cultures

Unité de vente :

24 unités
