

FICHE TECHNIQUE

Référence	354474
Désignation	Corning® BioCoat® Fibrillar Collagen



Description :

- ✓ Corning® BioCoat® Fibrillar Collagen Inserts with 1.0 μm PET Membrane in two 24-well Plates
- ✓ Brand Corning®
- ✓ Cell Growth Area 0,3 cm^2
- ✓ Membrane Material Polyethylene terephthalate (PET)
- ✓ Membrane Pore Size 3 μm
- ✓ Membrane Diameter 6,4 mm
- ✓ Packaging Format 12 Inserts/Plate/Pk, 2 Plates/Cs

Unité de vente :

24 unités
