

## FICHE TECHNIQUE

---

<b>Référence</b>	354515
<b>Désignation</b>	Plaque culture 6 puits



### Description :

- ✓ Corning® BioCoat™ Poly-L-Lysine 6-well Clear Multiwell Plate is comprised of clear TC-treated polystyrene with a uniform application of Poly-L-Lysine. Flat bottom
- ✓ Brand Corning®
- ✓ Cell Growth Area 9,6 cm<sup>2</sup>
- ✓ Storage 4° - 30°C
- ✓ Surface Coating Poly-L-Lysine

### Unité de vente :

5 unités

---