


## FICHE TECHNIQUE

---

<b>Référence</b>	354719
<b>Désignation</b>	Plaque 384 puits



---

### Description :

- ✓ Corning® PureCoat™ Amine 384-well Black Microplate
- ✓ No special storage or handling required
- ✓ Stable for at least 18 months from date of manufacture
- ✓ Quality control tested using an appropriate cell line
- ✓ Sterilized to SAL 10<sup>-6</sup> by gamma irradiation
- ✓ Noncytotoxic and nonpyrogenic
- ✓ Cultureware material is polystyrene suitable for cell culture
- ✓ Brand Corning®
- ✓ Storage Room Temperature
- ✓ Surface Coating Amine

### Unité de vente :

5 unités

---