

## FICHE TECHNIQUE

---

<b>Référence</b>	354773
<b>Désignation</b>	Plaques 6 puits

---

### Description :

- ✓ Store at room temperature — no special storage or handling required :
- ✓ Stable for at least 18 months from date of manufacture :
- ✓ Quality control tested using an appropriate cell line :
- ✓ Sterilized to SAL  $10^{-6}$  by gamma irradiation :
- ✓ Non-cytotoxic and nonpyrogenic :
- ✓ Cultureware material is polystyrene suitable for cell culture :
- ✓ Surface Coating : Carboxyl

### Unité de vente :

5 unités

---