

FICHE TECHNIQUE

Référence	354784
Désignation	Boite de petri



Description :

- ✓ Store at room temperature - No special storage or handling required
- ✓ Stable for at least 18 months from date of manufacture
- ✓ Quality control tested using an appropriate cell line
- ✓ Sterilized to SAL 10^{-6} by gamma irradiation
- ✓ Noncytotoxic and nonpyrogenic
- ✓ Cultureware material is polystyrene suitable for cell culture
- ✓ Brand : Corning®
- ✓ Cell Growth Area : 58,1 cm²
- ✓ Storage : Room Temperature
- ✓ Surface Coating : Carboxyl

Unité de vente :

10 unités
