

## FICHE TECHNIQUE

---

<b>Référence</b>	356241
<b>Désignation</b>	Corning® PureCoat™ Fibronectin Peptide

---

### Description :

- ✓ Corning® PureCoat™ Fibronectin Peptide 24-well Clear Plate, with Lid, Individually Wrapped, Sterile
- ✓ cGMP and 13485 compliant
- ✓ Brand Corning®
- ✓ Cell Growth Area 2 cm<sup>2</sup>
- ✓ Storage Room Temperature (15-30°C)
- ✓ Surface Coating Fibronectin
- ✓ Sterile Yes

### Unité de vente :

10 unités

---