

FICHE TECHNIQUE

Référence	356255
Désignation	Matrigel matrix pour culture d'organoïde 10ml

Description :

- ✓ Brand : Corning®
- ✓ Size : 10 mL
- ✓ Source : Engelbreth-Holm-Swarm lathrytic mouse tumor
- ✓ Storage Recommendation : Store at -20°C.
- ✓ Handling : Avoid multiple freeze-thaws. Do not store in frost-free freezer. KEEP FROZEN.
- ✓ Application : Air-Liquid Interface (ALI) Cultures
- ✓ Protein Concentration : Lot-dependent
- ✓ Cell Types : Organoid, Organoid Lines from Primary cells and Stem cells, Stem and progenitor cells
- ✓ Usage : Vial
- ✓ Surface Type : ECM
- ✓ Surface Category : Extracellular Matrices (ECMs) and Biological Coatings

Unité de vente :

1 unité
