

FICHE TECHNIQUE

Référence	356271
Désignation	Corning® PureCoat™ Collagen I

Description :

- ✓ Corning® PureCoat™ Collagen I Peptide 24-well Clear Plate, with Lid, Sterile
- ✓ cGMP and 13485 compliant
- ✓ Brand Corning®
- ✓ Cell Growth Area 2 cm²
- ✓ Storage Room Temperature (15-30°C)
- ✓ Surface Coating Collagen I
- ✓ Sterile Yes

Unité de vente :

10 unités
