

## FICHE TECHNIQUE

---

<b>Référence</b>	356273
<b>Désignation</b>	Flacon culture



---

### Description :

- ✓ Quality: cGMP and 13485 compliant
- ✓ 175cm<sup>2</sup> Rectangular Canted Neck Cell Culture Flask with Vented Cap
- ✓ Brand : Corning®
- ✓ Flask Shape : Rectangular
- ✓ Neck Style : Canted
- ✓ Cap Style : Vented
- ✓ Surface Coating : Collagen I

### Unité de vente :

10 unités

---