

CORNING	Life Sciences	Document Number: LSR00023	Rev.: 5
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets, 5 mL 0.10 mL Increments			Page: 1 of 2

APPLICATIONS:

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in an aseptic environment.

FEATURES:

- Guaranteed sterile claim (individual)
 - Contamination control
- Consistent performance
 - Reliability, accuracy is +/-2% at full volume
- Aseptic Presentation
 - Contamination control
- Non-pyrogenic
 - Low endotoxins
- Color coding
 - Ease of identification on packaging and blow-out rings
- Reverse graduations
 - Easy reading of the amount of liquid removed or the amount remaining
- Negative graduations
 - Can measure past the stated capacity
- Crisp, legible, rotogravure printing
 - Easily read sample volumes
- Crystal grade polystyrene
 - Can clearly see fluid sample
- Non-wettable
 - Minimum meniscus
- Individually wrapped packaging: "Thermoform" polyethylene-coated paper on one side, clear plastic on the other
 - Reduces "static cling", quick and easy to open. Two ways to open: Peel or Pop-through

STERILITY:

Product is gamma irradiated

RESTRICTIONS:

- Refer to "Thermoplastics Properties Chart," available at:
http://catalog2.corning.com/Lifesciences/media/pdf/an_DL_226_Falcon_ThermoPlastics_Broch.pdf

CORNING	Life Sciences	Document Number: LSR00023	Rev.: 5
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets, 5 mL 0.10 mL Increments			Page: 2 of 2

Falcon® 5 mL, 0.10 mL Increments, Serological Pipet

Material: Pipet - Polystyrene; Plug - Polyester fiber

Related Catalog Numbers:	Cat. No.	Packout Description
Sterile	356543	50/bag, 200/case
Sterile	357529	25/bag, 500/case
Sterile	357543	Individually wrapped, 200/case

Dimensions: = XXX mm (XXX inches)

TIP:

Length: = 37.36 (1.471)
 Interior Diameter: = 1.57(.062)
 Outside Diameter: = 3.35 (.132)

MOUTHPIECE:

Length: = N/A
 Interior Diameter: = 8.03 (.316)
 Outside Diameter: = 6.01 (.236)
 Plug Length: = 15.09 (.594)

EXTRUSION:

Interior Diameter: = 6.01 (.236)
 Outside Diameter: = 8.03 (.316)
 Length: = 311.15 (12.250)

PIPET SIZE: 5 mL

INCREMENTS: 1/10


GRADUATIONS: Reverse

OVERFILL CAPACITY: 2.5 mL

COLOR CODE: Blue

GENERAL DIMENSIONS:

A = 348.151(13.721)
 B = 3.35 (.132) Diameter
 C = 8.03 (.316) Diameter
 D = 8.03 (.316) Diameter



Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/lifesciences/trademarks.
 © 2013 Corning Incorporated