


FICHE TECHNIQUE

Référence	356650
Désignation	Corning® BioCoat® Collagen I



Description :

- ✓ Corning® BioCoat® Collagen I 96-well White/Opaque Flat Bottom TC-treated Microplate, with Lid
- ✓ Brand Corning®
- ✓ Plate Format 96-well
- ✓ Plate Color Clear bottom white
- ✓ Well Bottom Flat
- ✓ Well Bottom Color Clear
- ✓ Well Volume 392 µL
- ✓ Cell Growth Area 0,32 cm²
- ✓ Sterile No
- ✓ Lids Included Oui
- ✓ Surface Coating Collagen I

Unité de vente :

50 unités
