

## FICHE TECHNIQUE

---

<b>Référence</b>	356652
<b>Désignation</b>	Plaque culture 6 puits



---

### Description :

- ✓ Corning® BioCoat™ Gelatin 6-well Clear Multiwell Plate is comprised of clear TC-treated polystyrene with a uniform application of gelatin. Flat bottom
- ✓ Brand Corning®
- ✓ Cell Growth Area 9,6 cm<sup>2</sup>
- ✓ Storage 4° - 30°C
- ✓ Surface Coating Gelatin

### Unité de vente :

50 unités

---