

## FICHE TECHNIQUE

---

<b>Référence</b>	356664
<b>Désignation</b>	Plaque 384 puits



---

### Description :

- ✓ Corning® BioCoat™ Collagen I 384-well White/Clear Plate is comprised of white TC-treated polystyrene with a uniform application of rat tail collagen type I. Clear flat bottom, with lid, sterile
- ✓ Brand Corning®
- ✓ Storage 4° - 30°C
- ✓ Surface Coating Collagen I
- ✓ Sterile Yes

### Unité de vente :

50 unités

---