

FICHE TECHNIQUE

Référence	356719
Désignation	Plaque 384 puits



Description :

- ✓ Corning® PureCoat™ Amine 384-well Black Microplate
- ✓ No special storage or handling required
- ✓ Stable for at least 18 months from date of manufacture
- ✓ Quality control tested using an appropriate cell line
- ✓ Sterilized to SAL 10⁻⁶ by gamma irradiation
- ✓ Noncytotoxic and nonpyrogenic
- ✓ Cultureware material is polystyrene suitable for cell culture
- ✓ Brand Corning®
- ✓ Storage Room Temperature
- ✓ Surface Coating Amine

Unité de vente :

50 unités
