

FICHE TECHNIQUE

Référence	356773
Désignation	Plaques 6 puits



Description :

- ✓ Store at room temperature — no special storage or handling required
- ✓ Stable for at least 18 months from date of manufacture
- ✓ Quality control tested using an appropriate cell line
- ✓ Sterilized to SAL 10^{-6} by gamma irradiation
- ✓ Non-cytotoxic and nonpyrogenic
- ✓ Cultureware material is polystyrene suitable for cell culture
- ✓ Surface Coating : Carboxyl

Unité de vente :

50 unités
