	<b>Life Sciences</b>	Document Number: <b>LSR00021</b>	Rev.: <b>6</b>
Description (Class): <b>Reference, Customer Technical Data Sheet – Serological Pipets , 25 mL 0.25 mL Increments</b>		Page: <b>1 of 2</b>	

**APPLICATIONS:**

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in an aseptic environment.

**FEATURES:**

- Consistent performance
  - Reliability, accuracy is +/-2% at full volume
- Non-pyrogenic
  - Low endotoxins
- Color coding
  - Ease of identification on packaging and blow-out rings
- Reverse graduations
  - Easy reading of the amount of liquid removed or the amount remaining
- Negative graduations
  - Can measure past the stated capacity
- Crisp, legible, rotogravure printing
  - Easily read sample volumes
- Crystal grade polystyrene
  - Can clearly see fluid sample
- Individually wrapped packaging:
  - Two ways to open: Peel or Pop-through

**REGULATORY COMPLIANCE:**

Falcon® products are manufactured under the current ISO 9001 and ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

**STERILITY:**

Product is gamma irradiated

**RESTRICTIONS:**

- Refer to "Thermoplastics Properties Chart," available at:  
[http://catalog2.corning.com/Lifesciences/media/pdf/an\\_DL\\_226\\_Falcon\\_ThermoPlastics\\_Broch.pdf](http://catalog2.corning.com/Lifesciences/media/pdf/an_DL_226_Falcon_ThermoPlastics_Broch.pdf)

<b>CORNING</b>	<b>Life Sciences</b>	Document Number: <b>LSR00021</b>	Rev.: <b>6</b>
Description (Class): <b>Reference, Customer Technical Data Sheet – Serological Pipets , 25 mL 0.25 mL Increments</b>			Page: <b>2 of 2</b>

**Falcon® 25 mL, 0.25 mL Increments, Serological Pipet**

**Material:** Pipet - Polystyrene; Plug - Polyester fiber

**Related Catalog Numbers:**

Sterile  
Sterile  
Sterile

**Cat. No.**

356525  
357515  
357525

**Packout Description**

Individually wrapped, 50/bag, 200/case  
Individually wrapped, 20/bag, 200/case  
Individually wrapped, 50/bag, 200/case

**Dimensions: = XXX mm (XXX inches)**

**TIP:**

Length: = 39.24 (1.545)  
Interior Diameter: = 2.36 (.093)  
Outside Diameter: = 4.78 (.188)

**MOUTHPIECE:**

Length: = 25.91 (1.020)  
Interior Diameter: = 5.08 (.200)  
Outside Diameter: = 7.62 (.300)  
Plug Length = 14.30 (.563)

**EXTRUSION:**

Interior Diameter: = 13.43 (.529)  
Outside Diameter: = 15.16 (.597)  
Length: = 249.35 (9.817)

**PIPET SIZE:** 25 mL

**INCREMENTS:** 0.25


**GRADUATIONS:** Reverse

**OVERFILL CAPACITY:** 7.0 mL

**COLOR CODE:** Magenta

**GENERAL DIMENSIONS:**

A = 310.90 (12.240)  
B = 4.78 (.188) Diameter  
C = 7.62 (.300) Diameter  
D = 15.16 (.597) Diameter



Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks).  
© 2013 Corning Incorporated