

CORNING	Life Sciences	Document Number: LSR00021	Rev.: 6
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets , 25 mL 0.25 mL Increments			Page: 1 of 2

APPLICATIONS:

For tissue culture usage or any non-toxic solution that requires aliquoted or aspirating/dispensing in an aseptic environment.

FEATURES:

- Consistent performance
 - Reliability, accuracy is +/-2% at full volume
- Non-pyrogenic
 - Low endotoxins
- Color coding
 - Ease of identification on packaging and blow-out rings
- Reverse graduations
 - Easy reading of the amount of liquid removed or the amount remaining
- Negative graduations
 - Can measure past the stated capacity
- Crisp, legible, rotogravure printing
 - Easily read sample volumes
- Crystal grade polystyrene
 - Can clearly see fluid sample
- Individually wrapped packaging:
 - Two ways to open: Peel or Pop-through

REGULATORY COMPLIANCE:


Falcon® products are manufactured under the current ISO 9001 and ISO 13485 Standards, and the current FDA Quality System Regulation 21 CFR Section 820.

STERILITY:

Product is gamma irradiated

RESTRICTIONS:

- Refer to "Thermoplastics Properties Chart," available at:
http://catalog2.corning.com/Lifesciences/media/pdf/an_DL_226_Falcon_ThermoPlastics_Broch.pdf

	Life Sciences	Document Number:	Rev.:
		[{ و و و و و }]	6
Description (Class): Reference, Customer Technical Data Sheet – Serological Pipets , 25 mL 0.25 mL Increments			Page: 2 of 2

Falcon® 25 mL, 0.25 mL Increments, Serological Pipet

Material: Pipet - Polystyrene; Plug - Polyester fiber

Related Catalog Numbers:	Cat. No.	Packout Description
Sterile	356525	Individually wrapped, 50/bag, 200/case
Sterile	357515	Individually wrapped, 20/bag, 200/case
Sterile	357525	Individually wrapped, 50/bag, 200/case

Dimensions: = XXX mm (XXX inches)


TIP:
 Length: = 39.24 (1.545)
 Interior Diameter: = 2.36 (.093)
 Outside Diameter: = 4.78 (.188)

MOUTHPIECE:
 Length: = 25.91 (1.020)
 Interior Diameter: = 5.08 (.200)
 Outside Diameter: = 7.62 (.300)
 Plug Length = 14.30 (.563)

EXTRUSION:
 Interior Diameter: = 13.43 (.529)
 Outside Diameter: = 15.16 (.597)
 Length: = 249.35 (9.817)

PIPET SIZE: 25 mL
INCREMENTS: 0.25
GRADUATIONS: Reverse
OVERFILL CAPACITY: 7.0 mL
COLOR CODE: Magenta

GENERAL DIMENSIONS:
 A = 310.90 (12.240)
 B = 4.78 (.188) Diameter
 C = 7.62 (.300) Diameter
 D = 15.16 (.597) Diameter



Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures.

For a listing of trademarks, visit www.corning.com/lifesciences/trademarks.
 © 2013 Corning Incorporated