


TECHNICAL DATA SHEET

Part Number	362781
Description	Vacutainer CPT Mononuclear Cell Preparation Tube



Benefits :

- ✓ The BD CPT contains FICOLL and is a one-step system to provide a simple method for isolating Peripheral Blood Mononuclear Cells (PBMC) - lymphocytes and monocytes from whole blood. The tube provides a sample-to-yield solution, eliminating stages of blood transfer, preparation of FICOLL, and manual gradient separation.
- ✓ Waiting time is significantly reduced as centrifugation can be carried out with the brakes enabled.
- ✓ The tube can yield up to 15 million PBMCs with centrifugation times of 15 or 20 minutes.

Technical specifications :

- ✓ Material:
Tube: Soda lime type glass. Silicone inner wall
Cap: Conventional cap - Butyl elastomer - Silicone cap
- ✓ Dimensions : 13 x 100 mm
- ✓ Volume : 4 ml

Sales unit :

60 Vacutainer CPT Mononuclear Cell Preparation Tube