

Laboratory table with glass worktop


Description:

- ✓ H type frame manufactured of steel profiles 60x30x2mm and 40x40x3mm
- ✓ Powder painted
- ✓ Frame is assembled with bolts
- ✓ Height adjustable feet (20mm range)
- ✓ Worktop is made of tempered glass placed on melamine faced chipboard
- ✓ 2mm ABS edges all around the worktop
- ✓ Thickness of the worktop 31mm
- ✓ Color: light grey (similar to RAL 7035)
- ✓ NET weight 33 kg
- ✓ Max evenly distributed load on the table 200 kg

Dimensions :

- ✓ Length 600mm
- ✓ Depth 900mm
- ✓ Height 900mm
- ✓ Service space in the back of the table 150mm

Available other options and not standard modifications under request



Sample picture