

### Laboratory table with ceramic worktop

**Description:**

- ✓ H type frame manufactured of steel profiles 60x30x2mm and 40x40x3mm
- ✓ Powder painted
- ✓ Frame is assembled with bolts
- ✓ Height adjustable feet (20mm range)
- ✓ Worktop is made of large size ceramic tile placed on melamine faced chipboard
- ✓ 2mm ABS edges all around the worktop
- ✓ Thickness of the worktop 33mm
- ✓ Color: light grey (similar to RAL 7035)
- ✓ NET weight 27 kg
- ✓ Max evenly distributed load on the table 200 kg

**Dimensions :**

- ✓ Length 600mm
- ✓ Depth 600mm
- ✓ Height 900mm
- ✓ No service space in the back of the table

*Available other options and not standard modifications under request*



Sample picture