

Laboratory table with ceramic worktop

Description:

- ✓ H type frame manufactured of steel profiles 60x30x2mm and 40x40x3mm
- ✓ Powder painted
- ✓ Frame is assembled with bolts
- ✓ Height adjustable feet (20mm range)
- ✓ Worktop is made of large size ceramic tile placed on melamine faced chipboard
- ✓ 2mm ABS edges all around the worktop
- ✓ Thickness of the worktop 33mm
- ✓ Color: light grey (similar to RAL 7035)
- ✓ NET weight 55 kg
- ✓ Max evenly distributed load on the table 200 kg

Dimensions :

- ✓ Length 1200mm
- ✓ Depth 900mm
- ✓ Height 900mm
- ✓ Service space in the back of the table 150mm

Available other options and not standard modifications under request



Sample picture