

### Laboratory table with ceramic worktop

**Description:**

- ✓ H type frame manufactured of steel profiles 60x30x2mm and 40x40x3mm
- ✓ Powder painted
- ✓ Frame is assembled with bolts
- ✓ Height adjustable feet (20mm range)
- ✓ Worktop is made of 2 large size ceramic tiles placed on melamine faced chipboard
- ✓ 2mm ABS edges all around the worktop
- ✓ Thickness of the worktop 33mm
- ✓ Color: light grey (similar to RAL 7035)
- ✓ NET weight 65 kg
- ✓ Max evenly distributed load on the table 200 kg

**Dimensions :**

- ✓ Length 1500mm
- ✓ Depth 900mm
- ✓ Height 900mm
- ✓ Service space in the back of the table 150mm



Sample picture

*Available other options and not standard modifications under request*