



SPECIFICATION SHEET

Device Code : 37-400-03

Date : 15-09-2021

Decive Name : Scalpel Handle No. 3

A

Device Specification :

1. Drawing : See Attached Drawing
2. Size : 126 mm
3. Finish : Satin
4. Marking : Device shall be stamped 'CE' Pakistan and stainless.
Name / Logo if required by the Customer.
5. Intended Use : Device shall be used for replaceable blades for dissection.
6. Function/ Attributes : Hold blad properly.

B

Composition Specification :

1. Steel Type used : 410 Steel Type
2. Steel Weight used : 1 piece shall have a weight of 20 gm \pm 3.
3. Required Hardness : Nil
4. Furnace Temperature : N/A

C


Component / Specification :

Component Required : 1 Component. Material Type Used: 410 steel type.
Component Required : Material Type Used:

D

List of Applicable Standards:

- 1) ASTM Standard For Raw Material
- 2) ISO 13485:2016
- 3) ISO 9001:2015
- 4) 93/42/EEC Class I

Prepared by 
& Date 15.9.21

Approved by 
& Date 15/9/2021