

TECHNICAL DATA SHEET

Part Number	380586
Description	Transfer pipet



Benefits :

- ✓ Material : LDPE
- ✓ Standard bulb
- ✓ Molded-in graduation marks that will never fade or wash off
- ✓ Graduated to 2,25 ml / graduation increments 0,25 ml
- ✓ Packaged in a unique dispenser box

Certificate :

- ✓ Non sterile

Technical specifications :

- ✓ Capacity : 5 ml
- ✓ Drop : 21 ml
- ✓ Bulb draw 1,8 ml
- ✓ Dimensions (stem \varnothing x tip \varnothing x length) 7,2 x 3 x 150 mm

Sales unit :

500 Transfer pipets