

TECHNICAL DATA SHEET

Part Number	38061-1EA-HON
Description	Iodide Iodate concentrate

Technical specifications :

- ✓ iodide iodate concentrate, for 1l standard solution, 0,05 m i2 (0,1n)
- ✓ CAS N°7758-05-6
- ✓ no data available
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ GHS pictograms GHS08 ; GHS07
- ✓ H Statements H372 ; H319 ; H315
- ✓ P Statements P280 ; P260

Sales unit :

1 Each Iodide Iodate concentrate
