

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	38210-1EA-HON
<b>Description</b>	Sodium hydroxide concentrate

---

## Technical specifications :

- ✓ sodium hydroxide concentrate, for 1 l standard solution, 0,1 m naoh (0,1n)
- ✓ CAS N°1310-73-2
- ✓ Molecular weight 40.0 g/mol
- ✓ Molecular formula HNaO
- ✓ Density 1.07 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each
  
- ✓ GHS pictograms GHS05
- ✓ H Statements H290 ; H314
- ✓ P Statements P260 ; P280

## Sales unit :

1 Each Sodium hydroxide concentrate

---