

TECHNICAL DATA SHEET

Part Number	38215-1EA-HON
Description	Sodium hydroxide concentrate

Technical specifications :

- ✓ sodium hydroxide concentrate, for 1 l standard solution, 1,0 m naoh (1,0 n)
- ✓ CAS N°1310-73-2
- ✓ Molecular weight 40.0 g/mol
- ✓ Molecular formula HNaO
- ✓ Density 1.350 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ GHS pictograms GHS05
- ✓ H Statements H290 ; H314
- ✓ P Statements P260 ; P280

Sales unit :

1 Each Sodium hydroxide concentrate
