

TECHNICAL DATA SHEET

Part Number	38243-1EA-HON
Description	Sodium thiosulfate concentrate

Technical specifications :

- ✓ sodium thiosulfate concentrate, for 1l standard sol., 0,01 m $\text{Na}_2\text{S}_2\text{O}_3$ (0,01n)
- ✓ CAS N°7772-98-7
- ✓ Molecular weight 158.11 g/mol
- ✓ Molecular formula $\text{Na}_2\text{O}_3\text{S}_2$
- ✓ Density 1.0 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ P Statements P280

Sales unit :

1 Each Sodium thiosulfate concentrate
