

TECHNICAL DATA SHEET

Part Number	38250-1EA-HON
Description	Oxalic acid concentrate

Technical specifications :

- ✓ oxalic acid concentrate, for 1l standard solution, 0,05 m (cooh)₂ (0,1n)
- ✓ CAS N°144-62-7
- ✓ Molecular weight 90.03 g/mol
- ✓ Molecular formula C₂H₂O₄
- ✓ Density 1.02 g/cm³
- ✓ Flash point Not applicable
- ✓ Material form liquid

- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each

- ✓ GHS pictograms GHS05
- ✓ H Statements H318
- ✓ P Statements P280

Sales unit :

1 Each Oxalic acid concentrate
