

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	38280-1EA-HON
<b>Description</b>	Hydrochloric acid concentrate

---

## Technical specifications :

- ✓ hydrochloric acid concentrate, for 1 l standard solution, 0,1 m hcl (0,1n) 85ml
- ✓ CAS N°7647-01-0
- ✓ Molecular weight 36.5 g/mol
- ✓ Molecular formula HCl
- ✓ Density 1.030 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each
  
- ✓ GHS pictograms GHS05
- ✓ H Statements H290
- ✓ P Statements P234 ; P280

## Sales unit :

1 Each Hydrochloric acid concentrate

---