

# TECHNICAL DATA SHEET

---

<b>Part Number</b>	38294-1EA-HON
<b>Description</b>	Sulfuric acid concentrate

---

## Technical specifications :

- ✓ sulfuric acid concentrate, for 1 l standard solution, 0,5 m h<sub>2</sub>so<sub>4</sub> (1,0 n)
- ✓ CAS N°7664-93-9
- ✓ Molecular weight 98.08 g/mol
- ✓ Molecular formula H<sub>2</sub>O<sub>4</sub>S
- ✓ Density 1.180 g/cm<sup>3</sup> (20 °C)
- ✓ Flash point Not applicable
- ✓ Material form liquid
  
- ✓ Storage condition Room Temperature
- ✓ Package type Plastic Ampoule
- ✓ Package size 1 Each
  
- ✓ GHS pictograms GHS05
- ✓ H Statements H314 ; H290
- ✓ P Statements P260 ; P280

## Sales unit :

1 Each Sulfuric acid concentrate

---